ER System | CeraPost®

ER CeraPost® – root posts made of zirconium-oxide ceramics.

The type of ceramic material used for ER CeraPosts® has proven successful for many years in dental and medical applications due to its excellent corrosion behaviour, complete biocompatibility and outstanding mechanical loading strength. Compared to metal posts, the tooth coloured CeraPost® offers the advantage of an aesthetically perfect restoration. At the same time, its high mechanical loading strength allows CeraPost® to achieve a lasting and stable restoration without micro-leakage by using the adhesive technique.

Material data of CeraPost®

Zirconium oxide ceramics:


Composition:	94,9% ZrO ₂
	5,1% Y ₂ O ₃
Bending strength:	820 MPa
Modulus of elasticity:	200 GPa
Hardness (Knoop):	18.000 N/mm ²

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All teeth and remaining tooth structures that are restored or reinforced with CeraPost® can be used as prosthetic abutments and restored with ceramic crowns, partial crowns and veneers without impairing the natural tooth colour, thus achieving an aesthetically pleasing result.

Advantages of adhesive fixation:

Compared to the conventional method, the retention force can be nearly doubled when using CeraPost®. Why is that? The adhesive technique in combination with roughened root canals permits a chemo-mechanical union of the used materials without micro-leakage. This advantage becomes particularly apparent in pre-prosthetic stabilization, where even a small tooth structure can be sufficiently stabilized without significant loss of tooth substance.

[Fig.1] Diagram of an adhesive bond

Restoration of teeth
with partially destroyed crown [coronal degree of destruction 10 - 70%]

Restoration of coronally destroyed teeth
with a 2-piece ceramic build-up [coronal degree of destruction 70 - 100%]

Clinical sequence of the aesthetic restoration of an incisor:


1. Preoperative clinical status.
2. Preoperative radiological status.
3. Preparation of the post site and the retention box using the ER root canal reamer 196 and the root facer 120D and checking the proper fit of the post.
4. Impression taking for the indirect post/core build-up.
5. Next, the build-up is created in the laboratory. The restoration is adapted to the individual requirements of the tooth and divided in two parts, i.e. the core build-up and the root post. The prefabricated ceramic post is placed into the build-up.
6. Checking the proper fit of the build-up, starting with the individually adapted core build-up.
7. This is followed by checking the proper fit of the reconstruction together with the CeraPost®.
8. Mechanical conditioning of the canal using roughening the instrument 196D. Insertion of the ceramic post reconstruction by means of dual-hardening composite.
9. Final trimming of the post/core build-up and impression taking.
10. Post operative clinical status after integration of the definite allceramic crown.

● 231L12.050

● 439L12.070

● 232L12.090

● 233L12.110

Reference:

Butz F, Lennon Á, Heydecke G, Strub JR.: Survival rate and fracture strength of endodontically treated maxillary incisors with moderate defects restored with different post- and core systems: An in vitro study. Int J Prosthodont 14, 58-64 (2001).

Dérand P, Dérand T.: Bond strength of luting cements to zirconium oxide ceramics. Int J Prosthodont 13, 131-5 (2000).

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Libman WJ, Nicholls JJ.: Load fatigue of teeth restored with cast posts and cores and complete crowns. Int J Prosthodont 8, 155-61 (1995).

Morgano SM, Brackett SE.: Foundation restorations in fixed prosthodontics: Current knowledge and future needs. J Prosthet Dent 82, 643-57 (1999).

Nergiz I, Schmage P, Platzer U, McMullan-Vogel C.: Effect of different surface textures on retentive strength of tapered posts. J Prosthet Dent 78, 451-7 (1997).

Özcan M, Vallittu PK.: Effect of surface conditioning methods on the bond strength of luting cement to ceramics. Dent Mater 19, 725-31 (2005).

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Rosentritt M, Fürer C, Behr M, Lang R, Handel G.: Comparison of in vitro strength of metallic and tooth-coloured posts and cores. J Oral Rehabil 27, 595-601 (2000).

Sahafi A., Peutzfeld A., Asmussen E., Gottfredsen K.: Bond strength of resin cement to dentin and to surface-treated posts of titanium alloy, glass fiber, and zirconia. J Adhes Dent 5, 153-62 (2003).

Schmage P, Sohn J, Nergiz I, Özcan M.: Various conditioning methods for root canals influencing the tensile strength of titanium posts. J Oral Rehabil 31, 890-4 (2004).

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Standlee JP, Caputo AA.: Effect of surface design on retention of dowels cemented with a resin. J Prosthet Dent 70, 403-5 (1993).

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