
































Carte de formes
Fraises en carbure de tungstène
Ordering guide | Tungsten carbide



Vitesses maximales permises / Maximum permissible speed

■ =  100 t/min./rpm	● =  16.000 t/min./rpm	◆ =  90.000 t/min./rpm
□ =  500 t/min./rpm	▲ =  20.000 t/min./rpm	■ =  100.000 t/min./rpm
▣ =  800 t/min./rpm	△ =  25.000 t/min./rpm	▣ =  120.000 t/min./rpm
● =  1.000 t/min./rpm	◆ =  30.000 t/min./rpm	▣ =  140.000 t/min./rpm
○ =  5.000 t/min./rpm	◇ =  35.000 t/min./rpm	□ =  150.000 t/min./rpm
⊙ =  6.000 t/min./rpm	◇ =  40.000 t/min./rpm	▣ =  160.000 t/min./rpm
⊖ =  7.000 t/min./rpm	◇ =  50.000 t/min./rpm	+ =  300.000 t/min./rpm
⦿ =  8.000 t/min./rpm	◆ =  60.000 t/min./rpm	 450.000 t/min./rpm
● =  10.000 t/min./rpm	◇ =  70.000 t/min./rpm	sans indication
○ =  15.000 t/min./rpm	◆ =  80.000 t/min./rpm	if not indicated

Tous ces instruments appartiennent à :

Dispositif médical : Classe IIa
Organisme notifié : TÜV Rheinland (0197)
Fabricant : Gebr. Brasseler GmbH & Co. KG
03/2019 - 410337V2
Pour toute information complémentaire
se référer à la notice d'utilisation

Quelques-uns des produits ou désignations dans le texte sont protégés par la loi sur les droits d'auteur, les droits sur la marque ou par un brevet. L'absence éventuelle d'une référence spéciale ou le signe ® n'exclut pas l'existence d'une telle protection juridique.

Cette œuvre est protégée par la loi sur les droits d'auteur. Tous droits réservés, y compris ceux de la traduction, copie et de la réimpression partielle ou complète. Tous droits de reproduction, même d'extraits (photocopie, microfilm ou un autre procédé, ou la sauvegarde sur informatique de toutes sortes) sont réservés et assujettis à notre autorisation écrite.



Sauf modification de produit ou couleur ou de fautes d'impression.

Some of the products and designations mentioned in the text are trademarked, patented or copyrighted. The absence of a special reference or the sign ® should not be interpreted as the absence of legal protection.

This publication is copyrighted. All rights, also with regard to translation, reprint and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or processed using electronic systems in any form or by any means (photocopying, microfilm or other methods) without the written permission of the editor.

Colours and products subject to alterations. Printing errors excepted.

Fraises
Burs





































500 204 001003 ● H1SEM.204.010
500 204 001003 ● H1SEM.204.012
500 204 001003 ● H1SEM.204.014
500 204 001003 ● H1SEM.204.016
500 204 001003 ● H1SEM.204.018
500 204 001003 ● H1SEM.204.021
500 204 001003 ● H1SEM.204.023
500 204 001003 ● H1SEM.204.025
500 204 001003 ● H1SEM.204.027
500 204 001003 ● H1SEM.204.029
500 204 001003 ● H1SEM.204.031
500 205 001003 ● H1SEM.205.010
500 205 001003 ● H1SEM.205.012
500 205 001003 ● H1SEM.205.014
500 205 001003 ● H1SEM.205.016
500 205 001003 ● H1SEM.205.018
500 205 001003 ● H1SEM.205.021
500 205 001003 ● H1SEM.205.023
500 205 001003 ● H1SEM.205.027







500 204 001003 ● H1SE.204.008
500 204 001003 ● H1SE.204.010
500 204 001003 ● H1SE.204.012
500 204 001003 ● H1SE.204.014
500 204 001003 ● H1SE.204.016
500 204 001003 ● H1SE.204.018
500 204 001003 ● H1SE.204.021
500 204 001003 ● H1SE.204.023
500 204 001003 ● H1SE.204.025
500 204 001003 ● H1SE.204.027
500 204 001003 ● H1SE.204.029
500 204 001003 ● H1SE.204.031
500 205 001003 ● H1SE.205.010
500 205 001003 ● H1SE.205.014
500 205 001003 ● H1SE.205.018
500 205 001003 ● H1SE.205.023
500 205 001003 ● H1SE.205.027

500 314 001003 ● H1SM.314.014
500 204 001003 ● H1SM.204.010
500 204 001003 ● H1SM.204.012
500 204 001003 ● H1SM.204.014
500 204 001003 ● H1SM.204.016
500 204 001003 ● H1SM.204.018
500 204 001003 ● H1SM.204.021
500 204 001003 ● H1SM.204.023
500 205 001003 ● H1SM.205.010
500 205 001003 ● H1SM.205.014
500 205 001003 ● H1SM.205.018
500 205 001003 ● H1SM.205.023
500 314 001003 ● H1S.314.010
500 314 001003 ● H1S.314.012
500 314 001003 ● H1S.314.014
500 314 001003 ● H1S.314.016
500 314 001003 ● H1S.314.018
500 314 001003 ● H1S.314.021
500 314 001003 ● H1S.314.023
500 314 001003 ● H1S.314.027










500 204 001003 ● H1S.204.008
500 204 001003 ● H1S.204.010
500 204 001003 ● H1S.204.012
500 204 001003 ● H1S.204.014
500 204 001003 ● H1S.204.016
500 204 001003 ● H1S.204.018
500 204 001003 ● H1S.204.021
500 204 001003 ● H1S.204.023
500 204 001003 ● H1S.204.025
500 204 001003 ● H1S.204.027
500 204 001003 ● H1S.204.029
500 205 001003 ● H1S.205.010
500 205 001003 ● H1S.205.014
500 205 001003 ● H1S.205.018
500 205 001003 ● H1S.205.023
500 205 001003 ● H1S.205.027
500 314 001001 H1.314.005
500 314 001001 H1.314.006
500 314 001001 H1.314.007
500 314 001001 H1.314.008
500 314 001001 H1.314.009
500 314 001001 H1.314.010
500 314 001001 H1.314.012
500 314 001001 H1.314.014
500 314 001001 H1.314.016
500 314 001001 H1.314.018

500 314 001001 H1.314.021
500 314 001001 H1.314.023
500 314 001001 H1.314.027
500 316 001001 H1.316.010
500 316 001001 H1.316.012
500 316 001001 H1.316.014
500 316 001001 H1.316.016
500 316 001001 H1.316.018
500 316 001001 H1.316.023
500 204 001001 H1.204.005
500 204 001001 H1.204.006
500 204 001001 H1.204.007
500 204 001001 H1.204.008
500 204 001001 H1.204.009
500 204 001001 H1.204.010
500 204 001001 H1.204.012
500 204 001001 H1.204.014
500 204 001001 H1.204.016
500 204 001001 H1.204.018
500 204 001001 H1.204.021
500 204 001001 H1.204.023
500 204 001001 H1.204.027
500 205 001001 H1.205.010
500 205 001001 H1.205.014
500 205 001001 H1.205.016
500 205 001001 H1.205.018
500 205 001001 H1.205.023






Fraises à séparer les couronnes
Crown cutters




Fraise pour le retrait des amalgames
Amalgam remover



Fraises à finir
« Q »
Q-Finishers



H134Q.314.014



H375RQ.314.016



H379Q.314.018



H379Q.314.023 +



H246Q.314.009



H390Q.314.018 +



H50AQ.314.010



H135Q.314.014 +



H48LQ.314.012 +




H48XLQ.314.012 +


Fraises à finir
Finishing instruments



H41.314.014



H41.314.018



H41.314.023 +



H41.314.027 ■



H41.204.014 ■



H41.204.018 ■



H41.204.023 ■



H41.204.027 ■



H46.314.012




H46.314.014



H46.314.018



H47L.314.012




H47L.314.014



H281.314.009



H282.314.010



H282.314.012




H283.314.010 +




H283.314.012 +



H283.314.014 +



H283.204.012 ■



H283E.314.012 +




H284.314.014 +




H281K.314.012




H282K.314.014



H282K.314.016




H282K.204.014 ■




H282K.204.016 ■




H283K.314.016



H283K.314.021 +



H283K.204.016 ■




H283K.204.021 ■




H284K.314.018 +



H297.314.012 +



H336.314.016



H336.314.018



H375R.314.012 +




H375R.314.014 +




H375R.314.016



H375R.314.018



H375R.314.021 +




H375R.314.023 +




H375R.314.012 +






































H375R.314.018



H375R.314.021 +



H375R.314.023 +



**Retouche
du titane
en bouche**
*Intraoral working
on titanium*



● **H856G.314.016**



● **H856G.314.018**



● **H856G.314.020 +**



● **H847KRG.314.016**



● **H847KRG.314.018**



● **H847KRG.314.020 +**



● **H379G.314.023 +**

Endodontie
Endodontics



■ **H269QGK.314.012**



■ **H269QGK.314.016**



■ **H269QGK.315.016**



500 314 219295
■ **H269GK.314.012 +**



500 314 219295
■ **H269GK.314.016**



▲ **H1SML31.205.004**



▲ **H1SML31.205.006**




▲ **H1SML31.205.008**




▲ **H1SML31.205.010**




▲ **H1SML31.205.012**




▲ **H1SML31.205.014**




▲ **H1SML34.205.004**



▲ **H1SML34.205.006**




▲ **H1SML34.205.008**



▲ **H1SML34.205.010**



▲ **H1SML34.205.012**



▲ **H1SML34.205.014**



■ **EX1.310.007**



◇ **EX1S.204.007**



◇ **EX1.204.007**




◇ **EX1L.204.007**




■ **EX2.310.011**




◇ **EX2S.204.011**



◇ **EX2.204.011**



◇ **EX2L.204.011**



■ **EX3.310.013**



◇ **EX3S.204.013**



◇ **EX3.204.013**



◇ **EX3L.204.013**

**Chirurgie/
Implantologie**
*Surgery/
Implantology*




500 205 001291	H141.205.010	■
500 205 001291	H141.205.014	■
500 205 001291	H141.205.018	■
500 205 001291	H141.205.023	■
500 205 001291	H141.205.027	■
500 205 001291	H141.205.031	■
500 205 001291	H141.205.035	■
500 205 001291	H141.205.040	■
500 206 001291	H141.206.010	■
500 206 001291	H141.206.014	■
500 206 001291	H141.206.018	■
500 206 001291	H141.206.023	■
500 206 001291	H141.206.027	■
500 206 001291	H141.206.031	■
500 206 001291	H141.206.035	■
500 206 001291	H141.206.040	■
500 104 001291	H141.104.010	■
500 104 001291	H141.104.014	■
500 104 001291	H141.104.018	■
500 104 001291	H141.104.023	■
500 104 001291	H141.104.027	■
500 104 001291	H141.104.031	■

500 104 001291	H141.104.035	■
500 104 001291	H141.104.040	◇
500 104 001291	H141.104.045	◇
500 104 001291	H141.104.050	◇
500 105 001291	H141.105.010	■
500 105 001291	H141.105.014	■
500 105 001291	H141.105.018	■
500 105 001291	H141.105.023	■
500 105 001291	H141.105.027	■
500 105 001291	H141.105.031	■
500 105 001291	H141.105.035	■
500 105 001291	H141.105.040	◇
500 105 001291	H141.105.050	◇
500 104 001291	H141Z.104.010	■
500 104 001291	H141Z.104.014	■
500 104 001291	H141Z.104.018	■
500 104 001291	H141Z.104.023	■
500 104 001291	H141Z.104.027	■
500 104 001291	H141Z.104.031	■
500 104 001291	H141Z.104.035	■
500 104 001291	H141Z.104.040	◇
500 104 001291	H141Z.104.045	◇
500 104 001291	H141Z.104.050	◇

500 205 001298	H141A.205.014	■
500 205 001298	H141A.205.018	■
500 205 001298	H141A.205.023	■
500 205 001298	H141A.205.027	■
500 205 001298	H141A.205.031	■
500 205 001298	H141A.205.035	■
500 205 001298	H141A.205.040	◇
500 206 001298	H141A.206.010	■
500 206 001298	H141A.206.014	■
500 206 001298	H141A.206.018	■
500 206 001298	H141A.206.023	■
500 206 001298	H141A.206.027	■
500 206 001298	H141A.206.031	■
500 206 001298	H141A.206.035	■
500 206 001298	H141A.206.040	◇
500 104 001298	H141A.104.010	■
500 104 001298	H141A.104.014	■
500 104 001298	H141A.104.018	■
500 104 001298	H141A.104.023	■
500 104 001298	H141A.104.027	■
500 104 001298	H141A.104.031	■
500 104 001298	H141A.104.035	■
500 104 001298	H141A.104.040	◇
500 104 001298	H141A.104.050	◇

500 205 001298	H141AZ.205.014	■
500 205 001298	H141AZ.205.018	■
500 205 001298	H141AZ.205.023	■
500 205 001298	H141AZ.205.027	■
500 205 001298	H141AZ.205.031	■
500 205 001298	H141AZ.205.035	■
500 205 001298	H141AZ.205.040	◇
500 206 001298	H141AZ.104.010	■
500 206 001298	H141AZ.104.014	■
500 206 001298	H141AZ.104.018	■
500 206 001298	H141AZ.104.023	■
500 206 001298	H141AZ.104.027	■
500 206 001298	H141AZ.104.031	■
500 206 001298	H141AZ.104.035	■
500 206 001298	H141AZ.104.040	◇
500 206 001298	H141AZ.104.050	◇



500 314 415298	H254E.314.012	◇
500 314 415298	H254LE.314.012	◇
500 204 415298	H254E.204.012	◇
500 205 415298	H254E.205.012	◇
500 104 415298	H254E.104.012	◇
500 314 415296	H254.314.010	◇
500 314 415296	H254.314.012	◇
500 314 415296	H255E.314.012	◇
500 314 415296	H255E.316.012	◇
500 314 415296	H255E.204.012	◇
500 314 415296	H255E.104.012	◇
500 104 001298	H162SL.314.014	■
500 104 001298	H162SXL.314.014	■
500 104 001298	H162SL.204.014	◇
500 104 001298	H162SL.104.014	■



H163A.204.014 ■




H163A.104.014 ■



500 316 171007
H33L.316.010 ■



500 316 171007
H33L.316.012 ■



500 316 171007
H33L.316.016 ■



500 316 194007
H33R.316.016 ■



H162ST.314.016 ■



H162ST.204.016 ■



H162ST.104.016 ■



H162STZ.314.016 ■




H162STZ.204.016 ■



H162STZ.104.016 ■




500 314 408298
H162A.314.016 ■



500 204 408298
H162A.204.016 ■



500 205 408298
H162A.205.016 ■



500 104 408298
H162A.104.016 ■



H162AZ.314.016 ■



H162AZ.204.016 ■



H162AZ.205.016 ■




H162AZ.104.016 ■




500 314 408297
H162.314.016 ■



500 204 408297
H162.204.016 ■




500 205 408297
H162.205.016 ■




500 104 408297
H162.104.016 ■




H162Z.314.016 ■



H162Z.204.016 ■




H162Z.205.016 ■



H162Z.104.016 ■



500 314 210295
H267.314.016 ■



500 104 210295
H267.104.016 ■



500 314 199295
H269.314.016 ■




H269Q.314.016 ■



H269Q.315.016 ■



500 314 408295
H161.314.016 ■



500 104 408295
H161.104.016 ■




H166ST.104.021 ■




H166STZ.104.021 ■



500 204 409298
H166A.204.021 ■




500 205 409298
H166A.205.021 ■




500 104 409298
H166A.104.021 ■



H166AZ.204.021 ■




H166AZ.205.021 ■




H166AZ.104.021 ■



500 204 409297
H166.204.021 ■



500 205 409297
H166.205.021 ■



500 104 409297
H166.104.021 ■



H166Z.104.021 ■




500 104 410297
H167.104.023 ◆




500 316 150001
H207D.316.012 +

Péri-implantite
Periimplantitis



● H379.310.014 ◆



● H379.310.023 ◆




○ H379UF.310.014 ◆



○ H379UF.310.023 ◆



● H48L.310.014 ◆



● H48L.310.023 ◆



○ H48LUF.310.014 ◆



○ H48LUF.310.023 ◆

Komet France SA
18, rue Fourcroy · 75017 Paris
Tél. +33 (0) 1 43 48 89 90
Fax +33 (0) 1 43 48 90 65
info@komet.fr
www.komet.fr
