


PiezoLine

Scaler/Perio | LU8-LU11, LU63



— max. 2 mm



Indikation/Indication

PiezoLUX und Sonosoft/LUX

	Leistungsstufe/ Power Level 10	Leistungsstufe/ Power Level 7	Leistungsstufe/ Power Level 6	
Scaler	LU9 LU10 LU11	-	LU8	
Implantatprophylaxe Implant prophylaxis	-	1981.KA1 SF19821	-	
Paro/Perio	LU63 LU64 LU65	-	-	

Hinweis: hohe Durchflussrate (mind. 50 ml/min.)

Note: High flow rate (at least 50 ml/min.)